



HUGER
VIDEO ENDOSCOPY



VIDEO VETERINARY GASTROSCOPY

Type	AGVE-2100AL	AGVE-2100B	AGVE-2100P	GVE-2100	GVE-2100P
Optical System	Field of view 120°(direct)				
Depth of View	3-100mm				
Tip Deflection	180°(Up),180°(Down). 160°(Left),160°(Right),	180°(Up),90°(Down).100°(Left),100°(Right)			
Distal Rigid Diameter	Φ 12.9mm	Φ 9.8mm	Φ 9.0mm	Φ 9.8mm	Φ 9.0mm
Insertion Tube Diameter	Φ 12.9mm	Φ 9.6mm	Φ 9.0mm	Φ 9.6mm	Φ 9.0mm
Working Channel	3.7mm	2.8mm	2.2mm	2.8mm	2.2mm
Working Length	3.3M	1.5M		1.05M	

VIDEO CENTER

Type	VIS-2100	VIS-2100B
Lamp	21W Xenon Lamp	50W Xenon Lamp
Life of lamp	500 Hours	2500 Hours
Color Temperature	5900°K	
Pump	Both durable and quiet pump	
Function	White Balance, Enhancement, Adjusting of Color, Digital Zoom, Internal memory which can store to 4 image, Picture in picture	
Video Out	Composite video*1, Y/C video*2,	
LCD Monitor	10 inch Monitor	



HUGER Endoscopy Instruments Co.,Ltd.

No,1328,Yaobei Road, Jiuting, Songjiang District,201615,Shanghai China

Tel:+86-21-67691764 Fax:+86-21-67691721

E-mail: inform@huger.cn <http://www.huger.cn>