


Surgery Catalogue



For more information on our products and services -
EMAIL: sales@optiquemedical.co.za

NP - 6

Urethro-cystoscope set (Click)



1. Connecting to TIANSONG light guide cable
Connecting to OLYMPUS light guide cable
Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STRYKER light guide cable



Surgery Catalogue


For more information on our products and services -
EMAIL: sales@optiquemedical.co.za

NP - 7


Urethro-cystoscopy set (Lockable)



2



1. Connecting to TIANSONG light guide cable
Connecting to OLYMPUS light guide cable
Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STYKER light guide cable




Surgery Catalogue

For more information on our products and services -



EMAIL: sales@optiquemedical.co.za

NP - 2

Pediatric cystoscopy set (Lockable)



- 1.Connecting to TIANSONG light guide cable
Connecting to OLYMPUS light guide cable
Connecting to STORZ light guide cable
- 2.Connecting to WOLF light guide cable
- 3.Connecting to ACMI,STRYKER light guide cable



Surgery Catalogue


For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

NQ - 3

Urethrotomy set (Lockable)

4



A0011 12°
Endoscope
Ø4×302mm

1. Connecting to TIANSONG light guide cable
- Connecting to OLYMPUS light guide cable
- Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STRYKER light guide cable



A3102 Sheath with obturator 21Fr



A3101 Working element (passive)



A3103 Supplementary sheath



A3104 Lancet



A3105 Half-round knife



A3106 Hook shape knife



T1001 Adaptor



A2145 Sealing cap

T5010 Light guide cable
Ø4.5×2000mm




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

DQ - 4

Resectoscopy set (Lockable)



A0012 30°
Endoscope
ø4 × 302mm

1. Connecting to TIANSONG light guide cable
- Connecting to OLYMPUS light guide cable
- Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STRYKER light guide cable



A3165.3S Outer sheath 26.5Fr



A3163.3S
Inner Sheath with obturator 24.5Fr



A3160.4S
Working element (Passive)



T1001 Adaptor



M1015
High frequency cord (Monopolar)



1,3-DS-24
Loop 30°



4-DS-24
Roller



6-DS-24
Spike



3,1-DS-24
Spud 30°



M0012
H.F Generator (Multifunction)



T4002
Evacuator (Ellik)
T1001.6
Adaptor



T5010
Light guide cable
ø4.5 × 2000mm

Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za




DQ - 6

Resectoscopy set (Click)

6



1. Connecting to TIANSONG light guide cable
- Connecting to OLYMPUS light guide cable
- Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STRYKER light guide cable




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

DQ - 7

Resectoscopy set (Click)



1. Connecting to TIANSONG light guide cable
- Connecting to OLYMPUS light guide cable
- Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STRYKER light guide cable



Surgery Catalogue


For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

DQ - 8

Resectoscopy set (Click)

8



A0011 12°
Endoscope
ø4×302mm

1. Connecting to TIANSONG light guide cable
- Connecting to OLYMPUS light guide cable
- Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STRYKER light guide cable



A3165.2 Outer sheath 25.5Fr



A3163.2 Inner Sheath with obturator 22Fr



A3162 Working element (Passive)



T1001 Adaptor



M1001 High frequency cord (Monopolar)



A3111 Loop



A3122 Roller



A3123 Spike



A3124 Spud



M0012 H.F Generator (Multifunction)



T4002 Evacuator (Ellik)
T1001.5 Adaptor




T5010 Light guide cable
ø4.5×2000mm

Surgery Catalogue


For more information on our products and services -
EMAIL: sales@optiquemedical.co.za

DQ - 9

Pediatric Resectoscopy set (Lockable)



1. Connecting to TIANSONG light guide cable
- Connecting to OLYMPUS light guide cable
- Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STRYKER light guide cable



Surgery Catalogue

For more information on our products and services -
EMAIL: sales@optiquemedical.co.za

PS - 3


Lithotripsy set (Curved head)



10

PS - 4

Lithotripsy set (Straight)



PS - 3

PS - 4

Lithotripsy interchangeable parts




Surgery Catalogue

For more information on our products and services -
EMAIL: sales@optiquemedical.co.za


11

SS - 3

Ureterorenoscopy set



Ureterorenoscopy set options



Surgery Catalogue


For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

JP - 1

Percutaneous nephroscope set

12




Surgery Catalogue

For more information on our products and services -
EMAIL: sales@optiquemedical.co.za

Urology instruments options

13




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Urology instruments options



Male urethral dilator 8Fr-30Fr

A2113	8Fr	A2113.8	16Fr	A2113.16	24Fr
A2113.1	9Fr	A2113.9	17Fr	A2113.17	25Fr
A2113.2	10Fr	A2113.10	18Fr	A2113.18	26Fr
A2113.3	11Fr	A2113.11	19Fr	A2113.19	27Fr
A2113.4	12Fr	A2113.12	20Fr	A2113.20	28Fr
A2113.5	13Fr	A2113.13	21Fr	A2113.21	29Fr
A2113.6	14Fr	A2113.14	22Fr	A2113.22	30Fr
A2113.7	15Fr	A2113.15	23Fr		



Female urethral dilator 18Fr-42Fr

A2114.1	18Fr	A2114.4	27Fr	A2114.7	36Fr
A2114.2	21Fr	A2114.5	30Fr	A2114.8	39Fr
A2114.3	24Fr	A2114.6	33Fr	A2114.9	42Fr




A2116 F17-F40

Urethral filar sound dilator






A2131.2	Fistula needle 20Fr
A2131.3	Fistula needle 18Fr
A2131.4	Fistula needle 22Fr
A2131.5	Fistula needle 16Fr



A2139.1
A2139.2

Light suction aspirator tip 3x280mm
Light suction aspirator tip 5x280mm



A2101.1
Biopsy forceps, rigid
7Fr x 395mm



A2102.1
Foreign body forceps, rigid
7Fr x 395mm



A2103.1
Serrated forceps, rigid
7Fr x 395mm



T1001.1
Adaptor for model 6, 7 resectoscopy set
T1001.5
Adaptor for model 8 resectoscopy set



T1001.2
Adaptor for lithotrite



T1001.3
Adaptor for veress needle



T1001.4
Adaptor for WOLF cystoscope



T1001.6
Adaptor for STORZ resectoscopy set



T1009
Luer-lock adaptor



A3125
Electrode, needle



A3126
Electrode, ball



A3127
Electrode, straight loop
(double shaft monopolar)



A3129
Electrode, roller
(double shaft monopolar)



A3128
Electrode, loop
(double shaft monopolar)





A3130
Electrode, spud
(double shaft monopolar)

Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Laparoscopy instruments




Surgery Catalogue


For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Detachable operating forceps



Flexible and rotatable forceps



Surgery Catalogue

For more information on our products and services -


EMAIL: sales@optiquemedical.co.za

5 x 330 mm Laparoscopy instruments

17




E2005.1
Dissecting forceps, straight
(Head length 17mm)




E2005.2
Dissecting forceps, curved
(Head length 17mm)



E2005.3
Dissecting forceps, 30° angled
(Head length 17mm)



E2005.4
Dissecting forceps, 45° angled
(Head length 17mm)



E2005.5
Dissecting forceps, 70° angled
(Head length 17mm)



E2005.6
Dissecting forceps, 90° angled
(Head length 17mm)




E2005.7
Fixation grasping forceps, crocodile
(Head length 17mm)




E2005.8
Fixation grasping forceps, serrated
(Head length 13mm)




E2005.9
Fixation grasping forceps, loop
(Head length 17mm)




E2005.15
Grasping forceps



E2006.1
Scissors, straight
(Head length 13mm)



E2006.2
Scissors, straight, single action
(Head length 13mm)




E2006.3
Scissors, curved
(Head length 15mm)



E2006.4
Scissors, hooked, single action
(Head length 9mm)



E2006.10
Scissors, curved, head with ball
(Head length 16mm)



E2007
Grasping forceps
(Head length 14mm)



E2007.1
Blunt grasping forceps
(Head length 17mm)



E2007.4
Atraumatic grasping forceps
single action
(Head length 18mm)




E2007.5
Atraumatic grasping forceps
(Head length 15mm)




E2007.6
Grasping forceps,
rat toothed single action
(Head length 15mm)



E2007.7
Grasping forceps, double spoon
(Head length 10mm)



E2007.8
Gallbladder grasping
forceps, rough toothed
(Head length 18mm)



E2007.9
Gallbladder grasping forceps
(Head length 21mm)



E2007.10
Grasping forceps, greeve tooth
(Head length 10mm)



E2007.11
Grasping forceps
(Head length 21mm)



E2007.12
Grasping forceps, single action
(Head length 32mm)



E2007.27
Grasping forceps



E2007.33
Grasping forceps



E2007.34
Grasping forceps



E2007.35
Anastomosis traction forceps



E2007.36
Grasping forceps
(Head length 19mm)



E2007.39
Grasping forceps, rough tooth, hollow
(Head length 21mm)



E2007.45
Grasping forceps, reticulate



E2007.91
Gallbladder grasping forceps

Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

5 x 330 mm Laparoscopy instruments



E2007.13
Grasping forceps, flask style
(Head length 25mm)




E2007.14
Grasping forceps, clasp
(Head length 30mm)




E2007.16
Grasping forceps 2x3teeth, single action
with lock (Head length 30mm)




E2007.17
Grasping forceps 3x4teeth, single action
(Head length 30mm)



E2007.18
Atraumatic grasping forceps, Allie
(Head length 30mm)



E2007.19
Grasping forceps, Parker
(Head length 13mm)



E2007.20
Grasping forceps, single action
(Head length 30mm)




E2007.21
Grasping forceps
(Head length 22mm)




E2007.22
Appendix grasping forceps
(Head length 21mm)



E2007.23
Atraumatic grasping forceps
(Head length 30mm)



E2007.24
Grasping forceps, blunt-pointed with nail
(Head length 14mm)



E2007.25
Grasping forceps, heavy type
(Head length 18mm)



E2007.26
Grasping forceps 2x3teeth
(Head length 30mm)



E2007.31
Grasping forceps, round spoon head
(Head length 22mm)



E2007.32
Single tooth grasping forceps
(Head length 21mm)



E2018A6
Lung forceps
(Head length 32mm)




E2020 A
Tying forceps, type A
(Head length 17mm)




E2020.1 B
Tying forceps, type B
(Head length 17mm)




E2021
Small biopsy forceps, single action
(Head length 15mm)



E2022
Intestinal forceps
(Head length 22mm)




E2024
Grasping forceps, bullet head
(Head length 15mm)




E2024.1
Grasping forceps, round bullet head,
with spoon (Head length 15mm)



E2025
Grasping forceps, crocodile head
(Head length 15mm)



E2034
Grasping forceps, snake
(Head length 15mm)



E2035
Dissecting forceps
(Head length 30mm)



The code number for
flexible instruments
is ending with KW
Eg: E2005.1 - E2005.1KW



PPSU handle
(134° C 2bar)

1. Autoclave

2. The code number for PPSU handle series is beginning
with E9
Eg: E2005.1 - E9005.1

Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

10 x 330 mm Large laparoscopy instruments

19



E2002
Large grasping forceps, single action
(Head length 27mm)



E2003
Lithotomy forceps, single action
(Head length 30mm)



E2004
Large dissecting forceps, insulated
(Head length 22mm)



E2004.1
Large dissecting forceps, insulated
(Head length 22mm)




E2004.2
Large dissecting forceps, insulated
(Head length 22mm)



E2004.2A
Large dissecting forceps, insulated
(90°, head length 30mm)



E2004.4
Large dissecting forceps
(Head length 30mm)



E2004.6
Large dissecting forceps
(Head length 24mm)



E2004.5
Grasping forceps (Head length 60mm)



E2022.1
Intestinal forceps (Head length 68mm)



E2004.7
Grasping forceps, duck mouth
(Head length 30mm)



E2004.8
Grasping forceps, atraumatic
(Head length 36mm)



E2004.9
Large grasping forceps
(Head length 40mm)



E2005.410
Large dissecting forceps
(Head length 31mm)



E2004.11
Grasping forceps, atraumatic
(Head length 30mm)



E2004.12
90°
Large dissecting forceps



E2006
Large straight scissors
(Head length 18mm)



E2006.6
Large curved scissors
blade with teeth, single action
(Head length 24mm)



E2006.8
Large curved scissors
(Head length 28mm)



E2006.11
Large straight scissors



E2007.9C
Gallbladder grasping forceps
(Head length 35mm)



E2007.37
Grasping forceps, single tooth



E2018
Lung forceps
(Head length 43mm)




Surgery Catalogue


For more information on our products and services -

EMAIL: sales@optiquemedical.co.za


5 x 330 mm Bipolar coagulating forceps



5 x 330 mm Needle holder (head length 12mm)



3 x 330 / 250 / 220 mm Laparoscopy instruments




Surgery Catalogue


For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Electric coagulating instruments 5 x 330 mm



Suction and irrigation tube



Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Trocar set

22



With protection
E2016.3 ϕ 10.5
E2016.4 ϕ 6.5
E2016.7 ϕ 11



Spiral with protection
E2016.3A ϕ 10.5
E2016.4A ϕ 5.5



Cylindric
E2016.5 ϕ 10.5
E2016.6 ϕ 5.5
E2016.8 ϕ 11
E2016.9 ϕ 12
E2016.10 ϕ 6



Cylindric.spiral
E2016.5A ϕ 10.5
E2016.6A ϕ 5.5



Multifunctional valve
E2016.65 ϕ 10.5
E2016.66 ϕ 11



E2017.4 Reducer
(Engineering plastic) 10mm/5mm



E2017.5 Reducer
(Engineering plastic button) 10mm/5mm



Magnet
E2016.12 ϕ 5.5
E2016.13 ϕ 10.5
E2016.14 ϕ 11
E2016.15 ϕ 6
E2016.16 ϕ 12
E2016.19 ϕ 5.5

Stryker



Spiral magnet
E2016.17 ϕ 5.5
E2016.18 ϕ 10.5
E2016.24 ϕ 6



Magnet, with protection
E2016.20 ϕ 5.5
E2016.21 ϕ 6
E2016.22 ϕ 10.5
E2016.23 ϕ 11



Ball valve
E2016.41 ϕ 5.5
E2016.42 ϕ 6
E2016.43 ϕ 10.5
E2016.44 ϕ 11



Reducer (stainless steel)
E2017.1 10mm/5mm
E2017.2 5mm/3mm



E6004 Crossed sealing cap
5mm



E6005 Crossed sealing cap
10mm



E6001 Sealing cap 5mm




E6002 Sealing cap 10mm

Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Laparoscopy instruments




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

5 x 400 mm Single - hole umbilical laparoscopy instruments

24




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Single-hole umbilical laparoscopy instruments (curved)




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Mini laparoscopy instruments



Small incision instruments




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Thoracoscopy instruments

27




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Thoracoscopy instruments

28



Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

High frequency cord

29




Surgery Catalogue

For more information on our products and services -
EMAIL: sales@optiquemedical.co.za

ZCX

Automation suction and rinsing apparatus



30

Forceps



Ligation Forceps



Laparoscopy instruments, flexible head 10 x 330 mm



Surgery Catalogue

For more information on our products and services -
EMAIL: sales@optiquemedical.co.za

CJ - 1 Sigmoidoscopy set




CJ - 2 Rectoscopy set



CJ - 3 Sigmoidoscopy set



CJ - 1, CJ - 2, CJ - 3 Sigmoidoscopy & rectoscopy set interchangeable parts



Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

GQ - 6

Hysteroscopy set


32

F0107 0°
Hysteroscope $\phi 6.6 \times 220\text{mm}$
F0108 0°
Hysteroscope $\phi 5.3 \times 224\text{mm}$

F0101, 122°
Hysteroscope $\phi 6.6 \times 220\text{mm}$
F0106 22°
Hysteroscope $\phi 5.3 \times 220\text{mm}$



1. Connecting to TIANSONG light guide cable
- Connecting to OLYMPUS light guide cable
2. Connecting to STORZ light guide cable
2. Connecting to WOLF light guide cable
3. Connecting to ACMI, STYKER light guide cable



Flexible operating forceps $\phi 2.4 \times 395\text{mm}$



F2101
Biopsy forceps flexible




F2107
Loop catching forceps flexible



F2102
Scissors flexible



F2108
Scissors flexible



F2103
Grasping forceps with 4 claws flexible



Brush flexible
F2104 $\phi 2.4 \times 320\text{mm}$
F2104.1 $\phi 1.8 \times 320\text{mm}$



F2104 Brush rigid $\phi 2.5 \times 360\text{mm}$



F4001.2 Long needle



F4002
Fluid control uteromat (Automatic with computer)



A2145 Sealing cap



F4001
Fluid control uteromat (With infusion device and long needle)



F4001.1
Infusion device (Disposable)



T5010 Light guide cable
 $\phi 4.5 \times 2000\text{mm}$




Surgery Catalogue


For more information on our products and services -

EMAIL: sales@optiquemedical.co.za


GQ - 3 Hysteroscopy set



GQ - 5 Hysteroscopy set (Examine)



GQ - 7 Hysteroscopy set




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

GQ - 9


Hysteroscopy set (Examine)



34

GQ - 10

Hysteroscopy set (Examine)




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za


GQ - 11

Hysteroscopy set



GQ - 12

Hysteroscopy set




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

GQ - 8


Hysteroscopy set (Examine)



36

DQ - 5

Resectoscopy set




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Electric hysteria-cutter

37




Surgery Catalogue

For more information on our products and services -

EMAIL: sales@optiquemedical.co.za

Gynecology instruments

38



Surgery Catalogue

For more information on our products and services -
EMAIL: sales@optiquemedical.co.za

Cold light source



T5040 LED



T5001
PL-150 Halogen, double lamps 150W



T5002.1
GL-350E Xenon 350W



T5002.3
GL-50 Xenon 50W




T5004
PL-300 Halogen, double lamps



T5007.6
Portable cold light source (with charger)

Light guide cable diagram



Adaptor



T1003
Adaptor, for WOLF endoscope



T1004
Adaptor, for WOLF, TIANSONG cold light source



T1005
Adaptor, for OLYMPUS cold light source



T1006
Adaptor, for STRYKER light guide cable



T1007
Adaptor, for ACMI light guide cable



T1008
Adaptor, for ACMI, STRYKER endoscope

Surgery Catalogue

For more information on our products and services -


EMAIL: sales@optiquemedical.co.za

Endoscopy Camera System

NT-618

Endoscopy Camera System (WaterProof, LCD Panel)

- 1 618 Single chip medical camera (LCD Panel), 480lines, immersible camera lens from Japan
- 2 Coupler, Water Proof
- 3 15" SONY or 17" SAMSUNG Monitor
- 4 Pushcart, with multifunctional shelf



NT-501

Endoscopy Camera Documentation System

- 1 F-168D CCD camera, 480lines
- 2 Coupler, F23, Water Proof
- 3 Pushcart, with multifunctional shelf
- 4 Computer, Intel/I G/160G/DVD Recorder
- 5 Monitor, 19" LCD
- 6 Printer, HP4800DPI, high resolution print
- 7 Foot pedal
- 8 Video capture card, PHILIPS chip, highest quality of image, HD digital DVD level, dynamic magnification, dynamic filter processing
- 9 Software, can make DVD replay/record, print, many template/word stock/photo stock, powerful function of case management
10. Video teaching of the system, comprehensive teaching of the whole system
11. Automatic system recovery, fast recovery to factory setting in 20mins



NT-601

Endoscopy Camera Documentation System


- 1 F-068D CCD camera, 480lines
- 2 Coupler, F23, Water Proof
- 3 Pushcart, with multifunctional shelf
- 4 Computer, Intel/I G/160G/DVD Recorder
- 5 Monitor, 19" LCD
- 6 Printer, HP4800DPI, high resolution print
- 7 Foot pedal
- 8 Video capture card, PHILIPS chip, highest quality of image, HD digital DVD level, dynamic magnification, dynamic filter processing
- 9 Software, can make DVD replay/record, print, many template/word stock/photo stock, powerful function of case management
10. Video teaching of the system, comprehensive teaching of the whole system
11. Automatic system recovery, fast recovery to factory setting in 20mins



NT-668

NT-668 Endoscopy Camera System (WaterProof, HD)

- 1 668 HD medical camera, 1080P, immersible camera lens from Japan
- 2 Coupler, Water Proof
- 3 21" SONY medical HD monitor
- 4 Pushcart, with multifunctional shelf



NT-638

Endoscopy Camera System (Water Proof, LCD Touch Panel)

- 1 638 Three chip medical camera (LCD Touch Panel), 800lines, immersible camera lens from Japan
- 2 Coupler, Water Proof
- 3 21" SONY or 15" JVC Monitor
- 4 Pushcart, with multifunctional shelf

Surgery Catalogue

For more information on our products and services -

CALL: +27 (0) 10 222 0020 FAX: +27 (0) 86 562 7568 EMAIL: sales@hbl-distributors.co.za



Camera system



A2500d F-068D

Single chip, water proof

A2500f 618

Single chip, water proof, LCD panel

A2500e 638

Three chip, water proof, LCD touch panel

A2500h 668

waterproof, HD(1080P)LCD touch panel

Coupler

Coupler, Straight



Code	Focal length
A2513.1	F18
A2512.3	F22
A2513.3	F25
A2512.1	F27
A2517	F28
A2512.5	F32
A2512.6	F35

Medical instruments sterilizing box



A0000.1

For cystoscopy set
530X230X90mm



A0000.2

For resectoscopy set
420X275X90mm



B0000.3

For sinuscopy instruments
380X210X70mm



B0000.2

For sinuscopy set
360X255X90mm



F0000

For electric hystera-cutter
560X275X140mm



19000100

Endoscope. Instruments in common use
550X245X70mm



E0000

For laparoscope
450X90X45mm



B0000

For sinuscope
280X70X50mm



B0000.4

For sinuscope
255X71X36mm



B0000.5

For sinuscope, autoclavable
255X71X36mm



A0000.4

For cystoscopy, autoclavable
390X71X36mm